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Indifference and reliability guarantee unrivalled performance and economical service!

Easy, play-free and centric running of chucks. Handpiece with triple ball bearing. Grip of drills, mandrels and shafts increasing with load. All handpieces with standardized collet chucks.

Handpiece of high precision with high-speed ball bearing for Dental Laboratories and Industrial Use.

Chucking: Burs, mandrels and other fixing parts must run concentrically. Used or distorted shanks will damage the chuck. The fixing of the drills or grinding shafts is effected by fixing the ring-grip B and screwing up the ring-grip A, or vice versa.

Releasing: Drills are released with stationary motor by fixing ring-grip A (right hand) and simultaneously turning ring-grip A (left hand) one half to full turn clockwise, eventually by introducing a used shank in one of the two holes and unscrew the ring-grip A. Remove tools when you are encountering a resistance.

Maintenance: All ball bearings are live time lubricated so all other lubrication is not necessary. Only the chuck should sometimes be cleaned and lubricated.

Do not try to disassemble the handpiece, since no access to covered parts is necessary. When using high-speed motors an initial heat generation will automatically subside.

Do not connect the handpiece to running drive shafts.

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Dental GmbH