QUBE II

Original

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Schick GmbH

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Schick GmbH

D21518 QUBE II Original Instructions

Thank you for purchasing a high-quality technical device from SCHICK GmbH. We hope you will enjoy working with your new dental handpiece and wish you the best of success. We have prepared this operating manual in order to acquaint you with your new device and to provide the necessary information for operation and maintenance.

Project data:

Project designation: QUBE II
Project number: 9000_12_02_2015
Order:

Trade name: Schick QUBE II

Product name: QUBE II

Serial number: Axxx xxx

Machine number: 9xxx

Model: QUBE II SF, ST, SK

Type: QUBE II

Manufacturer: Schick GmbH

Authorized representative:

Schick GmbH

Lehenkreuzweg 12
D-88433 Schemmerhofen
Tel. +49 7356 95000
Fax +49 7356 950095
info@schick-dental.de
www.schick-dental.de

Lehenkreuzweg 12
D-88433 Schemmerhofen
Tel. +49 7356 95000
Fax +49 7356 950095
info@schick-dental.de
www.schick-dental.de

Established: 2017/07

Revision date: 2018/01

Contents

1	About this operating manual	6
1.1	General notes	6
1.2	Related documentation	6
1.3	Signs and symbols used	6
1.4	Structure of warnings	7
1.5	Technical terms and abbreviations used	8
2	General safety regulations	8
2.1	Principles	8
2.2	Intended use	8
2.3	Foreseeable misuse	8
2.4	Selection and qualification of the personnel	9
2.5	Workstations for operating personnel	9
2.6	Safety equipment	9
	Fire protection	9
2.7	Safety signs	10
	Explanation	10
2.8	Safety measures	10
	Personal protective equipment	10
2.9	Safety regulations	10
	General information	10
	During transport	10
	During installation	10
	During operation	11
	Care and control work	11
	Maintenance and repair work	11
	Structural modifications	11
2.10	Upgrading and conversion	11
3	List of contents	11

4	Technical description	12
4.1	Overview	12
	Function description	13
4.2	Interfaces.....	13
4.3	Characteristics of the tools	13
4.4	Name plate	14
5	Transport and storage	14
5.1	Transport.....	14
	Installation site requirements	14
5.2	Storage.....	14
	Storage site requirements.....	14
6	Installation and commissioning.....	15
6.1	Installation	15
	Required tools	16
	Safety precautions before installation	16
	Installation	16
	Installing the control unit on the suspension strip	16
	Fixing the power supply	17
6.2	Commissioning	18
	Establish power supply	18
7	Operating unit.....	19
7.1	Overview and position	19
8	Operating software.....	19
8.1	Operating concept and LCD display.....	19
8.2	Function overview	21
8.3	Basic operation of the control unit	21
8.4	Operating the handpiece and controller	22
8.5	QUBE II switch on	22
8.6	Tool change on handpiece	22
8.7	Operating modes.....	22

8.8	Desktop variant	22
8.9	Maximum speed limitation.....	23
8.10	Speed control function.....	23
	Activating cruise control:.....	23
8.11	Changing the handpiece rotation direction	24
8.12	Using the optional electric wax knife	24
	Switching on the electric wax knife:	24
	Changing the heat output:.....	24
	Switching off the electric wax knife:	24
8.13	QUBE II Switching off	25
9	Troubleshooting	25
9.1	How to proceed in case of malfunctions and errors	25
9.2	Fault and error messages	25
10	Care and control tasks for operators	26
10.1	General information.....	26
10.2	Care and inspection plan.....	26
10.3	Performing care and control work	27
	Dismantling and cleaning of the chuck	27
	Dismantling the motor cable	28
11	Technical data	29
12	Annex	30
12.1	Service address.....	30
12.2	Declaration of Conformity	32

1 About this operating manual

You must read this operating manual before using the **QUBE II** for the first time.

In particular, you must observe Section 2 "General safety regulations".

1.1 General notes

These instructions are intended to acquaint you with your **QUBE II** and provide information on its proper use.

It contains important tips on how to operate the **QUBE II** safely and correctly. Compliance will:

- Avoid dangers
- Keep repair costs and downtimes to a minimum
- Increase the reliability and service life of the product

This manual must be read and adhered to by each person tasked with working with the **QUBE II**.

In addition to this operating manual, the regulations on accident prevention and environmental protection applicable at the place of installation must also be observed.

1.2 Related documentation

The latest version of this manual as well as up-to-date product information is always available under www.schick-dental.de.

1.3 Signs and symbols used

The following signs and symbols are used in this manual:

- Activity symbol: The text after this sign describes handling instructions that must be carried out in the order indicated, from top to bottom.
- ✓ Result symbol: The text after this sign describes the result of an action.

Info symbol: Additional Information

1.4 Structure of warnings

Warning levels	Signal word	Use in case of ...	Possible consequences if the safety advice is not observed:
	DANGER	Personal injury (imminent danger)	Death or severe injuries!
	WARNING	Personal injury (potentially hazardous situation)	Death or severe injuries!
	CAUTION	Personal injury	Slight or minor injuries!

Tab. 1.1 Warning levels

The warnings are structured as follows:

- Pictogram with signal word corresponding to warning level
- Description of danger (type of hazard)
- Description of consequences of the danger (resulting hazards)
- Measures (activities) to prevent the danger

DANGER!

Type of danger (text)

Consequences of danger (text)

➤ Prevention of danger (text)

Warning signs Special safety instructions are provided at the relevant locations. These are marked with the following symbols.

General danger point

This symbol is indicated before activities that may result in personal injury or extensive property damage.

If the source of the danger is unambiguous, this symbol is preceded by one of the following symbols.

High voltage

This symbol is indicated before activities which pose a risk of electric shock, possibly with fatal consequences.

Hand injuries

This symbol is indicated before activities that pose a risk of hand injury.

1.5 Technical terms and abbreviations used

Abbreviation	Explanation
MH	Motor handpiece
SK	Knee control
SF	Foot control
ST	Table control

Tab. 1.2 Abbreviations used

2 General safety regulations

2.1 Principles

Use only with the maximum speeds specified by the manufacturer. Use only tools which are approved by the manufacturer/dealer for the application and performance data of **QUBE II**

Applications: Industry, commerce

Qualifications of personnel: Trainees and interns must be instructed by an experienced operator

The optional wax knife allows for thermal modelling of the waxes.

The device is operated manually, optionally via an operating unit (table control), knee plate, foot control, foot pedal

2.2 Intended use

QUBE systems are exclusively designed for universal use in dental laboratories for the dry processing of solid workpieces, such as crowns and bridges, as well as plastic and cast-models. Only tools that are intended for the processing of such materials may be used.

2.3 Foreseeable misuse

The system may not be used under the following circumstances:

- in potentially explosive atmospheres

- for medical use on the patient
- processing of wet materials
- for flammable and combustible materials
- any use other than that specified by the manufacturer is considered misuse
- do not clean the handpiece with compressed air

2.4 Selection and qualification of the personnel

Qualifications of personnel: Trainees and interns must be instructed by an experienced operator

2.5 Workstations for operating personnel

Dental workplace in dental laboratory/practice laboratory/dental practice.

2.6 Safety equipment

The **QUBE II** lab unit features an operating unit displaying the pre-set maximum speed directly in the field of vision of the user.

Fire protection

No specific requirements

2.7 Safety signs

Reference to rating plate on the control unit.

Explanation

Thoroughly read the operating manual before commissioning.

2.8 Safety measures

Do not wear long hair untied.

Always use extraction.

Personal protective equipment

Protective goggles

Extraction at the workplace.

2.9 Safety regulations

General information

Caution!

Risk of injury from sharp and/or rotating tools!

Wear safety goggles.

During transport

Transport or shipping only in the original packaging or a suitable shipping box.

During installation

The installation procedure may not impair access to the electricity supply of the power supply unit.

During operation

Operate only with a sufficiently dimensioned extraction system and while wearing proper protective clothing.

In left-hand rotation, the chuck may become loose when used with a standard handpiece!

Care and control work

- In accordance with the operating manual
- Trainees and interns must be instructed by an experienced operator

Maintenance and repair work

- Only by the manufacturer
- Only trained personnel with specialist knowledge of the respective device class may implement configuration and teach-in tasks.

Structural modifications

Structural changes to the product are not permitted.

2.10 Upgrading and conversion

A risk of burning on the hot model probes exists when using the optional electric wax knife.

3 List of contents

knee control	Art.no.	foot control	Art.no.
QUBE II knee controller	9516	QUBE II foot controller	9517
motorhandpiece with cable	9500	motorhandpiece with cable	9500
control satellite	9520	control satellite	9520
handpiece rack with two	9127	handpiece rack with two	9127
keys for changing chuck		keys for changing chuck	
power supply	9102	power supply	9102
velcro	9130	mains cable	2160
mains cable	2160	cable satellite	9124

suspension strip	9103	pulse emitter for suction	9060
screws – 2 pieces-	3170	velcro	9130
cable satellite	9124		
pulse emitter for	9060		
suction			

Table control	Art.No.	Optional	Art.No.
QUBE II table control	9518	foot switch on/off	6370/2
motorhandpiece with cable	9500	dynamic foot control	9440
control satellite	9520		
handpiece rack with two keys for changing chuck	9127		
power supply	9102		
velcro	9130		
mains cable	2160		
suspension strip	9103		
screws – 2 pieces-	3170		
cable satellite	9124		
pulse emitter for suction	9060		

4 Technical description

4.1 Overview

The QUBE II system is a modern workplace system designed for universal use in dental laboratories for the processing of crowns and bridges, as well as plastic and cast-models. The optionally available wax knife also allows for the modelling of waxes. QUBE II is available in three different control variants: knee control (SK), foot control (SF) and table control (ST).


Function description

The handpiece is used as a hand-held grinding or separating tool. The handpiece speed is limited through the operating unit and its touch-sensitive surface. Depending on the version, the speed is variably controlled via the knee plate, foot lever or the connected dynamic foot pedal (table control unit only).

4.2 Interfaces

The control units of the **QUBE II** feature the following interfaces to other machines:


Connection for controlling extraction systems. The pulse emitter for the extraction signal must be connected in order to operate the in **QUBE II** conjunction with an automatic workstation extraction system. This is plugged into the pulse transmitter for suction of the control unit as shown in the figure, and on the other end into the corresponding socket of the workplace extraction system. If the **QUBE II** is operated with a KaVo SMARTair Single Workstation Extraction unit or a Zubler AP 501 Extraction Point Opener, the optional control line Article no. 9229 can be used for connection and to establish the switching signal.

4.3 Characteristics of the tools

Only tools with perfect concentricity may be used. The shaft diameter may be 2.35 mm or 3 mm, depending on the chuck variant.

4.4 Name plate

1. Manufacturer
2. Type
3. Serial number
4. Power supply
5. Output
6. Mains frequency
7. Order/article no.
8. Pictogram: Disposal information
9. Pictogram: Observe the operating manual
10. CE marking

5 Transport and storage

5.1 Transport

If external damage to the packaging is visible upon delivery of the goods, this must be immediately communicated to the carrier and confirmed in writing. The damage must then be immediately reported to Schick GmbH.

Installation site requirements

Highly cooled appliances must be brought to room temperature before commissioning.
Danger of condensation.


5.2 Storage

Storage site requirements

In the original packaging, only indoors, protected from moisture

6 Installation and commissioning

6.1 Installation

1. control unit (foot, knee or table version)
- 2a. power supply
- 2b. mains cable
3. motorhandpiece with cable
- 4a. control satellite
- 4b. connection cable satellite
5. pulse emitter for suction

6. electric wax knife (optional)
7. dynamic foot control (optional for table version)
8. foot switch on/off (optional for table version)

Required tools


For knee control/table control: screwdriver, possibly a drill/cordless screwdriver for fastening the suspension strip.

Safety precautions before installation

Unplug the power cord before working on the electrical equipment.

Installation

Installation of the suspension strip for the knee control unit or hidden installation of the table control unit.

- 1. The suspension strip is provided with spacers, which serve as a stop for correct positioning of the knee control unit.
- To do so, position the suspension strip on the worktable as shown in the figure and fasten with the screws provided. After assembly, simply break off the spacers on the suspension strip.

Installing the control unit on the suspension strip


- 2. To mount the knee control device, place the recessed backside on the suspension strip and push it backwards until you can feel it latch in.

Fixing the power supply

Fasten the power supply with the velcro provided.

The included self-adhesive velcro can be used to fasten the power supply in the footwell of the worktable.

1.

2.


- 1. Remove one of the two protective foils from the velcro strip and adhere it to the power supply.
➤ 2. Remove the second protective foil from the velcro strip and adhere the power supply to a clean, dry and flat surface.

Avoid continuous tensile loads on the power supply caused by improperly laid cables!

6.2 Commissioning

Establish power supply

- Connect the motor handpiece to the "MI" socket. Connect the operating unit to the control unit with the cable.
- Optional: Insert the plug of the electrical wax knife into the control unit connector.
- Insert the AC adapter plug into the power socket of the control unit. Connect the power cord to the power supply.

Do not plug in the power cord until the handpiece and the operating unit have been correctly connected to the control unit (100 - 230 Volt, 50/60 Hz)

- ✓ The device switches on as soon as the power supply is established.

7 Operating unit


7.1 Overview and position

Location of the control and display elements

8 Operating software

8.1 Operating concept and LCD display


After switching on, the LCD display of the **QUBE II** operating unit informs you about all operating conditions of the device.

In particular, the following parameters are displayed:

- *Readiness indicator:* (Switched on/off) The display shows "OFF" when the **QUBE II** is switched off.
- *Ready for operation:* The display shows the preselected maximum speed and the direction of rotation.
- *Operating state:* The speed symbol is lit continuously when the motor is running and flashes when the motor is stopped.
- *Speed:* Indicates the pre-set maximum speed (x 1000)
- *Load/torque:* The bar display shows the power currently used by the motor when loaded. If all bars are shown and the handpiece is continually used, the safety switch will turn off the motor.
- *Electric wax knife:* An illuminated symbol indicates whether the optional electric wax knife is operational.
- *Electric wax knife adjustment:* Indicates the pre-set heating value
- *Speed control:* Shows the activated speed control
- *Direction of rotation:* Shows the current direction of rotation of the motor

8.2 Function overview

8.3 Basic operation of the control unit

The operating unit features five touch-sensitive keys in the areas marked by indentations, as well as a touch-sensitive rotary control (hereinafter referred to as "wheel"). To actuate, simply touch the touch-sensitive keys with your finger or, in the case of the wheel, guide the finger circularly within the recess.

Every actuation of one of the keys is displayed on the LCD display.

Please refer to Section 8.1 for more information about the individual functions.


8.4 Operating the handpiece and controller

8.5 QUBE II switch on

Briefly tap the symbol on the operating unit. Three short beeps confirm that the unit is switched on.

8.6 Tool change on handpiece

The chuck can be opened or closed by turning the handle on the handpiece.

Note:

You must switch off the motor before changing tools! In order to preserve the chuck's accuracy and service life, a tool or the protective pin provided from factory must always be clamped - even when the unit is not being operated.

Caution: To avoid bending of the tool shafts at high speeds, always insert the tools used as far as possible into the chuck in order to achieve the maximum holding force!

8.7 Operating modes

- Variable control of handpiece speed via knee plate, foot lever and/or dynamic foot control for the table control.
- Speed control (see also Point 9.6).

8.8 Desktop variant

If the QUBE II is operated without a foot pedal (refer to accessory page 12), the handpiece is started/stopped via the centre button inside the wheel (see Fig. 8.2)

- Set the desired speed via the wheel on the control satellite
- Start the handpiece via the O button.

You must observe the maximum permissible speed of the clamped tool before starting the handpiece!

8.9 Maximum speed limitation

If a speed is set via the wheel on the control satellite, this cannot be overridden with the operating elements (foot lever or knee plate). The speed is controlled dynamically with the control element, up to the speed set.


To use a speed range between 200 and 1000 rpm, the maximum speed must be limited to 1,000 rpm on the operating unit.

A flashing speed symbol  indicates that the motor is stopped. Once the motor is running, the speed symbol is lit continuously.

8.10 Speed control function


The speed control function allows you to operate the motor handpiece without having to continuously actuate the operating element (knee plate, foot lever or dynamic foot control).

Activating cruise control:

- ✓ Press the  key once. The  symbol appears on the display.
- ✓ Once the speed is steadily held via the control element for more than two seconds, the control element can be released and the handpiece will continue to run with the last used speed. To stop the handpiece, briefly actuate the control element once.

8.11 Changing the handpiece rotation direction

The QUBE II handpiece can be changed from right-hand to left-hand rotation via the control satellite. The default direction of rotation is right-hand.

Changing the direction of rotation:

- Briefly press the button. The symbol on the LCD shows the selected direction of rotation.


If the QUBE II handpiece is subjected to greater forces in left-hand rotation, the chuck may become loose. If the QUBE II is used by left-handers with left-handed tools, it is recommended you use special left-handed handpiece Art.No. 9002/08.

8.12 Using the optional electric wax knife

If the optional electric wax knife, Art.No. 9066, is connected to the control unit of the QUBE II it can also be controlled via the control satellite.

Switching on the electric wax knife:

- Briefly press the button. The LCD display switches to displaying the heat output.

- Use the wheel to set a heating value between 1 and 20 within three seconds.
- ✓ After three seconds, the display will change back to the normal speed, show the symbol and adjusted heating value.

Changing the heat output:

- Briefly press the button once and adjust the heating value using the wheel within three seconds.
- ✓ The wax knife switches off if the heating value is not changed.

Switching off the electric wax knife:

- Briefly press the button.
- ✓ The wax knife switches off after three seconds without further input.

8.13 QUBE II Switching off

Press the  symbol on the operating unit for more than two seconds. Three short beeps confirm that the unit is switched off.

9 Troubleshooting

9.1 How to proceed in case of malfunctions and errors

Error messages disappear immediately from the display as soon as the knee plate or foot lever is released or the fault is eliminated. This means there are no waiting times for re-starting.

Exception: If the table control is operated without a foot switch, all error messages must be reset by simultaneously pressing and holding the touch fields O +  + . The control unit can be completely reset by disconnecting from the power supply.

If a fault cannot be solved with the above procedure, please contact an authorised service partner or Schick directly.

9.2 Fault and error messages

The QUBE II control unit features intelligent control electronics that can detect potential faults and display them as fault codes on the control satellite. If a malfunction occurs, this is indicated by a red-lit display of the control satellite, a red flashing control unit and an error message in the format "Exyy".

Example of error message: E0804 error

E0804/ E0201	No handpiece connected	Check whether the handpiece is correctly connected - check plug with bayonet connector!
E0806	Overload cut-off function	Decrease contact pressure on workpiece.
E0802	Handpiece blocked	Check whether the chuck is closed. Check smooth running of the shaft.
	Broken cable	Replace the handpiece cable.
General error	The control unit is flashing red, no display on the control satellite, keys without function	Check connection cable from the control satellite, remove the mains plug and plug it back in again.
General error	Keys are not very responsive	Unplug the mains plug and plug it back in again. Do not touch the control satellite when restarting, so as to avoid miscalibration!

10 Care and control tasks for operators

10.1 General information

The QUBE II motorhandpiece is designed for maximum durability. Despite this, the chuck should be removed and cleaned from time to time, and the dirt underneath the handpiece tip removed with a brush.

Contamination on the plastic surfaces can be removed with a dry cloth.

Do not use aggressive cleaning agents!


10.2 Care and inspection plan

Care and inspection plan for QUBE II handpiece:

Interval	Care and inspection tasks	Remarks
Weekly	Dismantle the tip and remove dirt	Do not use compressed air!
Once a month	Dismantle the chuck, clean and lubricate	Lubricate with Art.No. 51/1

10.3 Performing care and control work


Dismantling and cleaning of the chuck

1. Remove the tools, unscrew the tip and remove the dust cap from the handpiece
2. Insert the chuck key (Art.No. 4115) into the open chuck, hold in place on the wrench grip surface of the shaft with the open-ended spanner (Art.No. 9188) and pull out the chuck with a sudden movement or by knocking on the chuck key, if necessary.
3. Remove the chuck from the shaft

- i** After thoroughly cleaning the chuck, slightly grease it on the outside with special grease (Art.No. 51/1) before re-installation.
- i** Tighten the chuck only slightly when re-installing! In normal operation (right-hand rotation), the chuck tightens itself.
- i** In order to preserve the chuck's accuracy and service life, a tool or the protective pin provided from factory must always be clamped - even when the unit is not being operated.
- i** Never clean the handpiece with compressed air!
- i** The corresponding chuck tools can be found on the underside of the handpiece rack.

Dismantling the motor cable

1. Unscrew the tip from the handpiece, remove the handle sleeve with ring and cable cap in the direction of the arrow
2. Unplug the connector of the motor cable from the handpiece

When installing the cable connector, make sure the plug-in contacts are correctly positioned!

11 Technical data

motorhandpiece:

speed range	200 – 60.000 rpm
torque max.	9 Ncm
diameter	29 mm
length	140 mm
weight without cable	180 g
weight with cable	320 g
Cooling system	self sealing system without cooling
drive system	brushless DC-Motor
concentricity	< 0,02 mm
chuck	Ø 2,35 und 3,0 mm (optional)
tool change	Quick release

control units:

Dimensions	Knee control	Foot control	Table control
Width:	95 mm	155 mm	155 mm
Height:	203 mm	75 mm	55 mm
Depth:	210 mm	220 mm	155 mm
Weight:	902 g	1.820 g	555 g
	Control satellite	Power supply	
Width:	84 mm	175 mm	
Heigth:	60 mm	35 mm	
Depth:	115 mm	70 mm	
Weight:	280 g	682 g	
Operating voltage:	100 – 240 Volt		
Output	160 Watt		

The oscillation total value during operation is below 2,5 m/s² (EN 28662).

12 Annex

12.1 Service address

Schick GmbH

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D-88433 Schemmerhofen

Tel.: +49 7356 9500-0

Fax: +49 7356 950095

E-mail: info@schick-dental.de

Internet: www.schick-dental.de

For service, please send your device directly to Schick GmbH at the above address.

Schick GmbH
Lehenkreuzweg 12
D-88433 Schemmerhofen
Telefon +49 7356 9500-0
Telefax +49 7356 9500-95
E-Mail info@schick-dental.de
Internet www.schick-dental.de

Schick GmbH reserves the right to change this documentation and the descriptions, dimensions and technical data contained herein without prior notice.

Schick GmbH assumes no responsibility for possible errors in this documentation. In no event shall Schick GmbH be liable for damages or consequential damages resulting from the use of this documentation or the hardware and software described herein.

The user is responsible for backing up and updating data. Schick GmbH accepts no liability for loss or manipulation of existing data.

We expressly point out that any kind of duplication - even for internal business purposes - is strictly prohibited. The contents may not be made available to third parties and cannot be used for other purposes.

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declare herewith, that the product

QUBE II consisting of
 QUBE II motorhandpiece 9500 in connection with
 QUBE II unit 9516, 9517, 9518 and
 QUBE II control satellite 9520 and
 foot switch (on/off) 6370/2
 dynamic foot control 9440

Is in conformity with the following provisions of Directive:

2006/42/EG	(machinery directive)
2014/30/EU	(EMC directive)
2011/65/EU	(RoHS)

Name and address of person in charge	Wolfgang Schick Lehenkreuzweg 12 D-88433 Schemmerhofen
---	--

Schemmerhofen, April 2017

W. Schick
director

Subject to technical modifications

This unit complies with the current VDE (German association of electrical technicians) regulations concerning safety and suppression.

These operating instructions should be readily accessible and are best kept close to the unit itself.

We would like to take this opportunity to advise you that a proper repair service and suitable qualified personnel are required for such highly developed technical equipment. SCHICK GmbH guarantees to carry out perfect repairs using original spare parts.

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© Copyright Schick GmbH
© Copyright di proprietà di Schick GmbH
© Los derechos de autor pertenecen a Schick GmbH
We, Schick GmbH