



Mobile

The flexible solution

Highest precision without grid connection

Powerful and mobile

The newly developed rechargeable handpiece Q Mobile from Schick combines two worlds for the first time: unrestricted precision and mobility owing to independence from power sockets.

Your benefits

- ▶ Highest precision and concentricity properties
- ▶ High performance of 5.5 Ncm and 30,000 rpm
- ▶ Quick-release clamping
- ▶ Battery capacity display
- ▶ Can be switched between left- and right hand motion
- ▶ Optional: connecting cable for separation of battery and handpiece

Q Mobile Art.-No. 9705

Scope of delivery:

Q Mobile motorhandpiece without cable,
control unit with integrated battery,
charging station, transport case



SCHICK
create the future

Mobile

Optional accessories:

- ▶ Waist bag (without content) Art.-No. 9654
- ▶ Connecting cable for Q Mobile handpiece Art.-No. 9604
- ▶ Connecting cable for Q Basic, Q Profi, QUBE LV Art.-No. 9604/1



Safe and practical -
Supplied in a transport case

Technical data

Number of revolutions	3.000-30.000 rpm
Torque	5,5 Ncm
Operating voltage	90-264 Volt
Operating frequency	50/60 Hz
Warranty	1 year
Overload protection	Electrical

Handpiece

Diameter	Max. 29 mm
Length without cable	120 mm
Weight	194 g
Cooling	Self sealing system without cooling
Drive system	Brushless DC-Motor
Concentricity	< 0,02 mm
Chucks	Ø 2,35 mm (standard) Ø 3,0 mm (optional)
Tool change	Quick release

Battery

Width	71 mm
Height	125 mm
Depth	43 mm
Weight	250 g
Capacity	3.500 mAh



Schick GmbH

Lehenkreuzweg 12
D-88433 Schemmerhofen
Telefon +49 7356 9500-0
Telefax +49 7356 9500-95
info@schick-dental.de
www.schick-dental.de