

Features of the compatible milling tools

S W1 is compatible with the following milling machines: Schick S3 Master, S2 Master, S2 Profi and S1 Basic

For milling machines that have a retrofitted milling spindle with the dimensions listed below, S W1 is also compatible.



Schick GmbH recommends converting the milling spindle to stainless steel shaft and stainless steel collet to prevent corrosion. Please contact the service of Schick GmbH at +49 7356 9500-0 or info@schick-dental.de

In order to be able to use the S W1 lighthead with light, when using the S1 Basic, the power supply lighthead Art.no. 2480/1 needed.