

QUBE long version




Ersatzteilliste für Motorhandstück Art.-Nr. 9002
List of spare parts for motorhandpiece art.-no. 9002



4125 Stift / pin




6642 Ø2,35mm Spannzange / chuck
(6643 Ø3,00mm) Spannzange / chuck
(6641 Ø1,60mm) Spannzange / chuck



6653/1 Spitze kpl. / tip cpl.



6652 Dichtkappe 2,35mm / sealing cap 2,35mm
(6652/1) Dichtkappe 3,00mm / sealing cap 3,00mm




6662/1 Griffhülse / grip sleeve




7225/2 Übergangsring / ring


Werkzeuge für Spannzangenwechsel
Tools to change chuck



6223 Gegenschlüssel /
counter stay wrench



4115 Spannzangenschlüssel /
chuck key



9225 Motorkabel / motor cable

